



## 1. Product Name

Fibre Optic Table

## 2. Product Code

22790

## 3. Colour

The plywood frame has a dark wood finish whilst the fibre optic light strands automatically cycle through various colours.

## 4. Brief Description

Large enough to keep necessities to hand with the added benefit of sensory activity. When not being used to hold drinks, reading glasses etc. on the top, the 100 fibre optic strands can be placed over the user's lap to provide soothing mood lighting that can be safely held, touched and explored.

## 5. Contents

- 1 x Fibre Optic Table (*the cup, coaster and telephone are for illustration purposes only*)

## 6. Snoezelen® Stimulations

- Sight

## 7. Best Use

Install indoors and either enjoy passively or with one of ROMPA®'s Wi Fi controllers for wireless interactivity. The Fibre Optic Table is especially effective when used with other items from the Multi-Purpose Reminiscence Furniture range (see below). This product is not intended to be sat or stood upon. This table will hold smaller items such as cups, magazines, plants, ornaments etc. but will not hold a heavy load.

## 8. Compatible ROMPA® Products

ROMPA® offers an extensive range of sensory and reminiscence products. Please visit our websites [www.rompa.com](http://www.rompa.com) and [www.olderadults.co.uk](http://www.olderadults.co.uk) for more information.

### **Multi-Purpose Reminiscence Furniture**

22793 Grandfather Glow Clock

22794 Multi-Purpose Television

22795 Illuminated Book Case

22792 Faux Flame 'Fireplace'

22789 Rocking Chair

22791 Infinity Table

20525 Nostalgia Music Centre



### **Wi Fi compatibility**

This unit is ROMPA® Wi Fi compatible, and therefore will respond wirelessly to the following remote controllers:

- 18870, 20402 - Talking Cubes
- 19847 – Deluxe 8 Colour Wirefree Controller
- 21006 – 8 Colour Wirefree Controller
- 19939 – Wi Fi Interactive Mat Switch Set
- 21917 – Wi Fi Microphone by ROMPA®
- 21918 – Wi Fi Override Transmitter
- 20848, 21098 – Sensory Magic®

## 9. Starting Up

- Choose a suitable location near to a wall socket.
- Plug in.

## 10. Detailed Description

Large enough to keep necessities to hand with the added benefit of sensory activity. When not being used to hold drinks, reading glasses etc. on the top, the 100 fibre optic strands can be placed over the user's lap to provide soothing mood lighting that can be safely held, touched and explored. The strands can be stroked, whilst their weight over the knees is in itself calming. The strands change colour softly in sequence to maintain interest. With an integral shelf for additional storage of books, magazines and so on. Mains-powered.

Weight: 12kg

Size: 40L x 50W x 40cmH; 100cmL fibre optic harness

The Light Source used in this product turns on automatically once plugged into the mains, automatically cycling through a rainbow of colours. The Light Source is compatible with ROMPA® Wi Fi controllers and will change colour when a Wi Fi controller is in operation. When controlled in this way, the Light Source will change to the selected colour for one minute, before returning to automatic mode.

## 11. Safety

- **Always supervise the use of this product**
- This product weighs approximately 12kg – take care when lifting the table.
- This product is decorative and not intended to be sat or stood upon. This product is unlikely to withstand the weight of a person.
- For indoor use only, in a dry and clean environment.
- Make sure the power lead and the fibre optic strands are located safely out of the way and that these do not create a trip hazard.
- The fibre optic strands can be handled with care, but should not be handled roughly, pulled, mouthed, chewed, bent or trodden upon.

## 12. Technical Specification

Size: 40L x 50W x 40cmH  
Power: Input: 100-240v AC 50/60Hz 0.25A  
Output: 12v DC 1.25A 15W max  
CE: Meets the requirements of BS EN 9001-2000  
UK plug: 13 amp  
Weight: approximately 12kg

## 13. Installation

1. Choose an appropriate location for the Fibre Optic Table, bearing in mind that:
  - It needs to be plugged in (one plug)
  - The lighting effect will be strongest in darker environments than in bright day light
  - The product weighs approximately 12kg – please choose a suitably flat and stable floor
  - This product is free-standing.
2. Plug the product in. There is an ON/OFF switch on the Light Source. Once switched on, the fibre optic strands will light up. The light will change colour in sequence unless a ROMPA® wireless controller is used to choose a particular colour of light. The automatic colour light sequence will resume one minute after the ROMPA® wireless controller has been activated.
3. Suitably small, safe and light objects can be placed on the table if desired, as well as on the lower integral shelf.

## 14. Care and Maintenance

1. Periodically dust the product with a soft and dry cloth.
2. To remove any dirt or more persistent marks, use a soft and very slightly dampened cloth. Avoid the use of any cleaning products, especially abrasive cleaning products, bleach or anti-bacterial sprays as these can damage varnished/lacquered surfaces.
3. Treat this table as you would any other coffee table. For example, use coasters or mats when placing cups, mugs or warm crockery on the table.

## 15. Troubleshooting

**Problem:** *The fibre optic strands do not light up*

- Check that the power lead is plugged in
- Check that the Light Source is switched on
- Check that the fibre optic boss is still secured inside the Light Source circular port – reattach the boss of the fibre optic spray to the Light Source as required using the screw in the port

Further copies of these instructions can be downloaded at [www.rompa.com](http://www.rompa.com)

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at

[producthelp@rompa.com](mailto:producthelp@rompa.com)